

Contents

About the Company.....	1
Intensive Care Ventilator MV300	2
Intensive Care Ventilator MV200	4
High Frequency Jet Ventilator JV100 B.....	6
Depth of Anesthesia and Sedation Monitor MGA-06.....	7
Patient Monitors MPR6-03 for Intensive Care.....	8
Patient Monitors MPR6-03 for Anesthesiology.....	10
Patient Monitors MPR6-03 for Neonatology.....	12



Triton Electronic Systems Ltd.
is a unique Russian enterprise of full cycle:

- R&D Activities;
- Intensive Care and Anesthesia Equipment Manufacturing;
- OEM Medical Solutions.

All the processes on the enterprise are certified in accordance with international standard of Quality Management System ISO 13485.

Decades ago a large scientific and technical cluster was formed in the Urals. The region of Russia has accumulated scientific knowledge and R&D expertise in different sectors – from light industry to aerospace. Whole research institutes were relocated to the Urals, so new ideas and inventions amounted to thousands. As a result, the Ural Scientific School was founded. It trained excellent engineers and researchers capable of solving technical and scientific problems, develop high-precision instruments.

The company was founded in 1989 by biomedical researchers, engineers and physicians from Ural Scientific School. Our greatest achievements are closely connected with the name of MD Prof. Boris Zislin, who was the Chief anesthesiologist of Ural region for 37 years. Due to his cutting edge ideas we develop innovative medical technologies and manufacture a range of medical devices for anesthesiology and intensive care.

We are able to bring in unique and principally new solutions for global healthcare market. They are easier, more convenient and more reliable; moreover, they are the solutions that the market needs.




Boris D. Zislin, Professor, Doctor of Medicine, Honored doctor of Russian Federation.

One of founders of Ural anesthesiology and resuscitation school.

His research interest was high frequency jet ventilation (HFJV), he had made every effort to implement HFJV into the practice.

B. Zislin is the author of over 200 articles and four monographies.



Intensive Care Ventilator **MV300**

2



Intended use:

Ventilator is intended for controlled and assisted artificial ventilation of lungs for all patient groups with tidal volume from 1 ml in resuscitation units, surgery and intensive care departments, and also at transportation within professional medical facilities.

Categories of patients: adults, children, newborn infants, including pre-matures.

Display: 15", touchscreen, color, viewing angle adjustment.

Gas supply: built-in turbine; oxygen from central gas pipeline, cylinder or O₂ concentrator.

Power: 100–250 V, 50/60 Hz, built-in battery, up to 6 h of operation.

Interface: USB, Ethernet.

Ventilation parameters

Tidal volume, ml	1–3000
Breathing rate, bpm	1–120
T _{insp} , sec	0.2–10
Flow trigger sensitivity, lpm	0.1–20
Pressure trigger sensitivity, cmH ₂ O	0.1–20
PEEP, cmH ₂ O	0–50
Inspiratory pressure, cmH ₂ O	0–100
Pressure support, cmH ₂ O	0–80
I:E ratio	1:99–60:1

Ventilation modes

Mandatory ventilation	volume controlled	CMV VCV
	pressure controlled	CMV PCV
	pressure controlled, volume guaranteed	PCV VG
Synchronized intermittent mandatory ventilation	volume controlled	SIMV VC (+PS)
	pressure controlled	SIMV PC (+PS)
	double-controlled	SIMV DC (+PS)
Spontaneous breathing	with continuous positive airway pressure	CPAP (+PS)
	two levels of PEEP	BiSTEP (+PS)
	airway pressure release ventilation	APRV (+PS)
Noninvasive ventilation	noninvasive ventilation with a mask	NIV
	nasal CPAP	nCPAP (+PS)
	high flow oxygen therapy	HF_O ₂
Adaptive ventilation	intelligent support ventilation with auto-control of ventilation	iSV
Reserve mode	apnea / backup ventilation	Apnea / backup
Combined ventilation	combined high-frequency jet and conventional ventilation	+HFJV

Features and advantages

Patient age groups: adults, children, newborn infants, including prematures	Can be applied to all groups of patients: newborn infants, children and adults. The lowest tidal volume for neonates is 1 ml with digital proximal flow sensor. Built-in adult/pediatric flow sensor does not require replacement during the whole shelf life
Neonatal option	Specially designed ventilation mode for the extremely low weight neonates. Tidal volume is from 1 ml with ultrasensitive proximal flow sensor (from 0.1 lpm). The digital proximal flow sensor is very reliable and has a long life time. nCPAP ventilation mode with variable flow is available for the ventilation of neonates
Combined mechanical ventilation	Consolidates the advantages of jet and conventional mechanical ventilation. Jet ventilation enhances oxygenation of arterial blood in complicated situations and convective mechanical ventilation provides normal removal of carbon dioxide
AutoETS	Auto-detection of the optimal expiratory trigger sensitivity and time of spontaneous inspiration in support ventilation modes (CPAP, NIV, BiSTEP) provides patient comfort and reduces the work of breathing. It can be switched off and the desired sensitivity level of the expiratory trigger can be set manually
AutoRAMP	Auto-detection of pressure rise time allows to match the ventilator flow to the patient's respiratory needs
Evaluation of metabolic needs	The method of indirect calorimetry. Continuous measurement: oxygen consumption ($\dot{V}O_2$), carbon dioxide elimination ($\dot{V}CO_2$), respiratory quotient (RQ), resting energy expenditure (REE)
Advanced monitoring of alveolar ventilation parameters	Volumetric capnography ($\dot{V}CO_2$), airway dead space volume (V_d), alveolar minute ventilation volume (\dot{V}_{alv})
Auxiliary pressure monitoring	Measurement of esophageal pressure/trachea via esophageal catheter and transpulmonary pressure measuring channel via balloon catheter
Stress index	Integral parameter indicating correctness of PEEP and tidal volume VT settings
Adaptive ventilation (iSV)	Mode of intelligent adaptive ventilation supporting patients with any level of breathing activity. It automatically determines the parameters of control and support ventilation pressure. This mode is ideal from weaning of patient from the ventilator
Mainstream capnometry	Monitoring of $EtCO_2$, $FiCO_2$ and respiratory rate. Quick response QuRe™ sensor provides the highest precision of capnogram waveform for true clinical diagnostics and allows working even in HF ventilation mode up to 200 bpm
Integrated calculator of cardiac output by Fick method	Calculation of the cardiac output parameter (CO) based on the data of alveolar ventilation monitoring
Pulse oximetry	Measurement of oxygen saturation of arterial blood hemoglobin SpO_2 , pulse rate (PR), photoplethysmogram

Intensive Care Ventilator **MV200**



Intended use:

Ventilator is intended for controlled and assisted artificial ventilation of lungs for all patient groups with tidal volume from 10 ml in resuscitation units, surgery and intensive care departments, and also at transportation within professional medical facilities.

Categories of patients: adults, children, infants.

Display: 12.1" or 15", touchscreen, color, viewing angle adjustment.

Gas supply: built-in turbine; oxygen from central gas pipeline, cylinder or O₂ concentrator.

Power: 100–250 V, 50/60 Hz, built-in battery, up to 6 h of operation.

Interface: USB, Ethernet.


Ventilation parameters

Tidal volume, ml	10–3000
Breathing rate, bpm	1–120
T _{insp} , sec	0.2–10
Flow trigger sensitivity, lpm	0.5–20
Pressure trigger sensitivity, cmH ₂ O	0.5–20
PEEP, cmH ₂ O	0–50
Flow _{max} [†] , lpm	180
Inspiratory pressure, cmH ₂ O	0–100
Pressure support, cmH ₂ O	0–80
I:E ratio	1:99–60:1

Ventilation modes

Mandatory ventilation	volume controlled	CMV VCV
	pressure controlled	CMV PCV
	pressure controlled, volume guaranteed	PCV VG
Synchronized intermittent mandatory ventilation	volume controlled	SIMV VC (+PS)
	pressure controlled	SIMV PC (+PS)
	double-controlled	SIMV DC (+PS)
Spontaneous breathing	with continuous positive airway pressure	CPAP (+PS)
	two levels of PEEP	BiSTEP (+PS)
	airway pressure release ventilation	APRV (+PS)
Noninvasive ventilation	noninvasive ventilation with a mask	NIV
Adaptive ventilation	intellectual support ventilation with auto-control of ventilation	iSV
Reserve mode	apnea / backup ventilation	Apnea / backup

Monitoring of ventilation parameters



Graphical monitoring:

- Curves: simultaneous displaying up to three curves (flow, pressure, volume, CO₂, SpO₂, VCO₂, P_{aux}, iSV).
- Monitoring parameters: simultaneous displaying up to 8 windows depending on user's choice.
- Loops: pressure–volume, flow–volume, flow–pressure, volume–auxiliary pressure.
- Controls: 8 parameters are displayed depending on ventilation mode.
- Hot keys: 9 icons for the fast access to special functions of ventilator.

Basic monitoring:

- Peak inspiratory pressure, plateau pressure, average pressure, PEEP, auto PEEP.
- Minute volume of respiration.
- Inspiratory, expiratory volume.
- Respiratory rate, spontaneous breaths frequency.
- Compliance C.
- Resistance R.
- I:E ratio.
- Concentration of inspired oxygen FiO₂.
- Leakage.
- Maximum flow on inspiration.
- EtCO₂, FiCO₂.

Advanced monitoring:

- AutoPEEP, total PEEP.
- Inspiratory time constant, expiratory time constant.
- Stress index.
- Index of respiratory effort.
- Work the patient's breathing, work of ventilator breathing.
- Inspiratory time.
- Coefficient of spontaneous breathing.
- Resistance to exhalation.
- Dynamic compliance.
- Circuit resistance.
- Circuit compliance.
- Shallow breathing index.
- SpO₂, pulse rate.
- Cardiac output.
- Volume of alveolar ventilation ($V_{alv\ min}$).

Features and advantages

Evaluation of metabolic needs	The method of indirect calorimetry. Continuous measurement: oxygen consumption (VO ₂), carbon dioxide elimination (VCO ₂), respiratory quotient (RQ), resting energy expenditure (REE)
Advanced monitoring of alveolar ventilation parameters	Volumetric capnography (VCO ₂), airway dead space volume (V _d), alveolar minute ventilation volume (MV _{alv})
Auxiliary pressure monitoring	Measurement of esophageal pressure/trachea via esophageal catheter and transpulmonary pressure measuring channel via balloon catheter
Stress index	Integral parameter indicating correctness of PEEP and tidal volume VT settings
Adaptive ventilation (iSV)	Mode of intellectual adaptive ventilation supporting patients with any level of breathing activity. It automatically determines the parameters of control and support ventilation pressure. This mode is ideal from weaning of patient from the ventilator
Mainstream capnometry	Monitoring of EtCO ₂ , FiCO ₂ and respiratory rate. Quick response QuRe™ sensor provides the highest precision of capnogram waveform for true clinical diagnostics and allows working even in HF ventilation mode up to 200 bpm
Integrated calculator of cardiac output by Fick method	Calculation of the cardiac output parameter (CO) based on the data of alveolar ventilation monitoring
Pulse oximetry	Measurement of oxygen saturation of arterial blood hemoglobin SpO ₂ , pulse rate (PR), photoplethysmogram

High Frequency Jet Ventilator **JV100 B**

6

Intended use: high-frequency jet ventilation with extended monitoring of ventilation and respiratory parameters.

Patient types: adults, children with body weight from 15 kg.

Display: 12.1" LCD, touchscreen, color, viewing angle adjustment.
Gas mixer and FiO₂ control.

Mainstream CO₂ sensor.

Power: 220 V/50 Hz, built-in battery, from 1 h of operation.

Gas supply: oxygen from central gas pipeline, cylinder; air – from central gas pipeline or compressor.

Humidification: build-in automatic heating system with thermal control and humidification.

Ventilation parameters

FiO ₂ control in non-injection and catheter modes, %	21–100
Respiratory rate, bpm	30–300
Minute volume (modes w/ or w/o Venturi jet injector), lpm	5–30
Minute volume (catheter mode), lpm	3–20
I:E ratio	1:3; 1:2; 1:1; 3:2; 2:1

Ventilation modes

High frequency jet	injection	HFJV
	catheter	
	non-injection	

Additional functions:

- The pulmonary percussion.
- The expiratory pause.
- The oxygenation.

Monitoring of artificial pulmonary ventilation parameters:

- Monitoring of minute volume.
- Monitoring of respiratory rate.
- Monitoring of the carbon dioxide concentration.
- Monitoring of peak inspiratory pressure.
- Monitoring of positive end expiratory pressure.
- Ratio monitoring of inhalation duration to exhalation duration I:E.
- Monitoring of carbon dioxide partial pressure.
- Monitoring of oxygen concentration.
- Monitoring of the minimum, maximum and average airway pressure.
- Monitoring of the inspiration volume.

Graphic display:

- Flow graphic.
- Volume graphic.
- Pressure.
- Capnogram.



Depth of Anesthesia and Sedation Monitor

MGA-06



Specially designed for patient depth of anesthesia and sedation monitoring.

7

Technical specification

Display 5", touchscreen, color

Built-in battery (2 h of operation)

24 hours of graphical trends

Visual and sound alarm system

Universal mounting on a horizontal surface, in vertical or horizontal poles

Power supply: 100–250 V, 50/60 Hz

Interface: Wi-Fi

Main parameters

AI – Brain activity index

indicates patient's depth of anesthesia and sedation level. Depth of anesthesia evaluation is based on analysis of electroencephalogram (EEG). Algorithm of EEG analysis keeps information about typical signs of influence of different groups of anesthetics – types of EEG suppression

SR – EEG signal suppression rate

reflects the suppression of EEG activity and determines level of the further deepening of anesthesia

SQI – EEG signal quality index

indicates the influence of impedance of cable electrodes, artifacts, high-frequency noise, power mains disturbances and decreases linearly from 100% to 0%

EMG – Electromyographic component level

electrical activity of facial muscles

Features and advantages

Universal

Operates with inhalation and intravenous anesthetics

Low cost in operation

Kit from 3 standard ECG-electrodes is used for measurement

Easy to use and safe for patient

3 ECG-electrodes are placed on the fronto-temporal lobes of the patient

Efficient

Using the device allows:

- reduce the consumption of anesthetics, as well reduce the awakening from anesthesia time by 35–50%;
- stably maintain the required depth of sedation (according to statistics, without AI-control more than 69% of patients have insufficient or excessive depth of sedation);
- mostly eliminate the risk of anesthesia awareness.



90–100	awake
80–89	light sedation
60–79	coma, sedation
40–59	surgical stage
20–39	deep anesthesia, burst-suppression patterns
0–19	extremely deep anesthesia

Convenient

Monitor is light and portable. It can be placed on the working surface or fixed on the rail or rack

Patient Monitors MPR6-03

Intended use: Intensive Care

8



Designed to monitor the main parameters of patient's vital signs.

Applications:

- Patient age groups: adults, children and infants.
- In operating rooms and intensive care units.
- In the wards of anesthesiology departments.
- In the specialized departments (neurology, cardiology and others).
- In the departments and offices of functional diagnostics.
- In the general admission departments of hospitals.

Vital signs and modules

Electrocardiography (ECG)	12 leads: I, II, III, aVL, aVR, aVF, V ₁ -V ₆ Heart rate variability analysis (HRV) ST-analysis Heart rate Analysis and automatic detection of arrhythmia
Pulse oximetry TREATON®	Pulse rate SpO ₂ Photoplethysmogram
Pulse oximetry Masimo SET® 	Pulse rate SpO ₂ Photoplethysmogram
Thermometry	Two channels Temperature difference measurement
Impedance method of measuring respiration parameters	Respiratory rate Respirogram
NIBP	Sys / Mean / Dia arterial pressure
Automatic continuous measurement of blood pressure (cNIBP)	Sys / Mean / Dia arterial pressure, pulse rate transit time relative to QRS
Noninvasive monitoring of central hemodynamics parameters	CO, SV, CI, SVR, SVRI, SVI, LCW, LCWI, LSW, LSWI, LSP, LSPI, VSV Original technology of volume-compression oscillometry
Invasive measurement of blood pressure (2 IBP channels)	Sys / Mean / Dia blood pressure Selection of pressure type: ART, PA, CVP, ICP, RAP, LAP, RVP, UA IBP1, IBP2 waveforms
Minimally invasive measurement of CVP and other low blood pressures in various areas of the body (IIND 500/75)	Measurement range, mmH ₂ O: -200...450
Mainstream capnography	EtCO ₂ , FiCO ₂ , RR Capnogram
Sidestream capnography	EtCO ₂ , FiCO ₂ , RR Capnogram
Sidestream oximetry	EtO ₂ , FiO ₂
Continuous monitoring of metabolic parameters (based on data from the respiratory mechanics module)	Automatic calculation of oxygen consumption and CO ₂ elimination with respiratory mechanics VCO ₂ , VO ₂ , REE, RQ
Respiratory mechanics	PIP, PEEP, Ve, MVe, flow and pressure waveforms Volumetric capnogram
Metabolic calculator	VO ₂ , VCO ₂ , REE, RQ
Depth of anesthesia (sedation level)	Brain activity index AI Electroencephalogram (EEG)

**Expert technologies:**

- Continuous monitoring of the metabolic needs of the patient.
- Respiratory mechanics.
- Noninvasive monitoring of central hemodynamic parameters.
- Capnography in the mainstream.

Basic functions

Touchscreen	12.1" or 15"
Power supply	220 V/50 Hz
Battery	2 h of operation
Trends	digital and graphic, 240 h
Interface	wired or wireless connection to the central monitoring station (up to 32 monitors) wired connection to the hospital network USB for transferring patient data to PC
Alarm system	three levels of priority setting alarm limits alarm history
Data printing	three channels thermal printer
Medical calculators	for drug doses for oxygenation and ventilation parameters for renal function for hemodynamics parameters

Benefits

Metabolic monitoring	<p>Emergency and planned surgery is always accompanied by the development of hypermetabolism-hypercatabolism syndrome. Late and unadequate nutritional support of the critically ill patient after surgery leads to the development of a number of adverse effects and complications. Principles of early clinical nutrition and maximum adequacy to the specific metabolic needs of the operated patient are a good decision for successful healing process</p> <p>The method of indirect calorimetry without use of disposables. Easy to use – method based on gas analysis data</p> <p>Continuous measurement:</p> <ul style="list-style-type: none"> • oxygen consumption (VO_2); • carbon dioxide elimination (VCO_2); • respiratory quotient (RQ); • resting energy expenditure (REE).
Noninvasive monitoring of central hemodynamic parameters	<p>Cardiac output assessment by a non-invasive volumetric compressive oscillometry (VCO). VCO bases on measurements of the blood volume change under an inflatable cuff around a limb. VCO determines blood pressure by recording the volumetric arterial oscillograms allowing speed and character of pulse-wave increase and decrease evaluation, and judgment about heart pump condition and vascular walls elasticity</p>

Patient Monitors MPR6-03

Intended use: Anesthesiology

10




Designed to monitor the main parameters of patient's vital signs.

Applications:

- Patient age groups: adults, children and infants.
- In operating rooms and intensive care units during anesthesia care.
- In the wards of anesthesiology departments.
- In the specialized departments (neurology, cardiology and others).

Vital signs and modules

Electrocardiography (ECG)	12 leads: I, II, III, aVL, aVR, aVF, V ₁ -V ₆ Heart rate variability analysis (HRV) ST-analysis Heart rate Analysis and automatic detection of arrhythmia
Pulse oximetry TREATON®	Pulse rate SpO ₂ Photoplethysmogram
Pulse oximetry Masimo SET® 	Pulse rate SpO ₂ Photoplethysmogram
Thermometry	Two channels Temperature difference measurement
Impedance method of measuring respiration parameters	Respiratory rate Respirogram
NIBP	Sys / Mean / Dia arterial pressure
Automatic continuous measurement of blood pressure (cNIBP)	Sys / Mean / Dia arterial pressure, pulse rate transit time relative to QRS
Noninvasive monitoring of central hemodynamics parameters	CO, SV, CI, SVR, SVRI, SVI, LCW, LCWI, LSW, LSWI, LSP, LSPI, VSV Original technology of volume-compression oscillometry
Invasive measurement of blood pressure (2 IBP channels)	Sys / Mean / Dia blood pressure Selection of pressure type: ART, PA, CVP, ICP, RAP, LAP, RVP, UA IBP1, IBP2 waveforms
Minimally invasive measurement of CVP and other low blood pressures in various areas of the body (IIND 500/75)	Measurement range, mmH ₂ O: -200...450
Mainstream capnography	EtCO ₂ , FiCO ₂ , RR Capnogram
Continuous monitoring of metabolic parameters (based on data from the respiratory mechanics module)	Automatic calculation of oxygen consumption and CO ₂ elimination with respiratory mechanics: VCO ₂ , VO ₂ , REE, RQ
Respiratory mechanics	PIP, PEEP, Ve, MVe, flow and pressure waveforms Volumetric capnogram
Metabolic calculator	VO ₂ , VCO ₂ , REE, RQ
Anesthetic agent module (multigas)	Fi/Et: CO ₂ , O ₂ , N ₂ O, 5AA. RR, MAC Selection of waveform: CO ₂ , O ₂ , anesthetic agents
Depth of anesthesia (sedation level)	Brain activity index AI Electroencephalogram (EEG)



Expert technologies:

- Anesthetic agent module (multigas).
- Depth of anesthesia (sedation level).
- Respiratory mechanics.
- Noninvasive monitoring of central hemodynamic parameters.
- Capnography in the mainstream.

Basic functions

Touchscreen	12.1" or 15"
Power supply	220 V/50 Hz
Battery	2 h of operation
Trends	digital and graphic, 240 h
Interface	wired or wireless connection to the central monitoring station (up to 32 monitors) wired connection to the hospital network USB for transferring patient data to PC
Alarm system	three levels of priority setting alarm limits alarm history
Data printing	three channels thermal printer
Medical calculators	for drug doses for oxygenation and ventilation parameters for renal function for hemodynamics parameters

Benefits

Depth of anesthesia and sedation level monitoring

Analyzing of patient's EEG in time and frequency domains and defines brain activity index (AI), which reflects patient's depth of anesthesia and sedation. Based on analysis of electroencephalogram (EEG). Algorithm of EEG analysis keeps information about typical signs of influence of different groups of anesthetics – types of EEG suppression



- 90–100 awake
- 80–89 light sedation
- 60–79 coma, sedation
- 40–59 surgical stage
- 20–39 deep anesthesia, burst-suppression patterns
- 0–19 extremely deep anesthesia

Anesthesia agent analysis – continuous sidestream measurement and monitoring of gas concentration in the patient's airways

Measured parameters:

- concentration of CO₂ at inspiration (FiCO₂) and at the end of expiration (EtCO₂);
- concentration of O₂ at inspiration (FiO₂) and at the end of expiration (EtO₂);
- concentration of N₂O at inspiration (FiN₂O) and at the end of expiration (EtN₂O);
- concentration of anesthetic gas agents (Ax): Halothane (Hal), Enflurane (Enf), Isoflurane (Iso), Desflurane (Des), Sevoflurane (Sev) at inspiration (FiAx) and expiration (EtAx);
- respiration rate (RR).

Patient Monitors MPR6-03

Intended use: Neonatology

12




Designed to monitor the main parameters of patient's vital signs.

Applications:

- Monitoring of newborn infants with extremely low body weight from 500 g.
- In neonatal operating rooms and intensive care units.
- In pediatric departments (neurology, cardiology and others).
- In admission pediatric departments.

Vital signs and modules

Electrocardiography (ECG)	12 leads: I, II, III, aVL, aVR, aVF, V ₁ -V ₆ Heart rate variability analysis (HRV) ST-analysis Heart rate Analysis and automatic detection of arrhythmia
Pulse oximetry TREATON®	Pulse rate SpO ₂ Photoplethysmogram
Pulse oximetry Masimo SET® 	Pulse rate SpO ₂ Photoplethysmogram
Thermometry	Two channels Temperature difference measurement
Impedance method of measuring respiration parameters	Respiratory rate Respirogram
NIBP	Sys / Mean / Dia arterial pressure
Automatic continuous measurement of blood pressure (cNIBP)	Sys / Mean / Dia arterial pressure, pulse rate transit time relative to QRS
Invasive measurement of blood pressure (2 IBP channels)	Sys / Mean / Dia blood pressure Selection of pressure type: ART, PA, CVP, ICP, RAP, LAP, RVP, UA IBP1, IBP2 waveforms
Minimally invasive measurement of CVP and other low blood pressures in various areas of the body (IIND 500/75)	Measurement range, mmH ₂ O: -200...450
Mainstream capnography	EtCO ₂ , FiCO ₂ , RR Capnogram
Sidestream capnography (sample rate 50 ml/min)	EtCO ₂ , FiCO ₂ , RR Capnogram
Sidestream oximetry	EtO ₂ , FiO ₂
Metabolic calculator	VO ₂ , VCO ₂ , REE, RQ
Anesthetic agent module (multigas)	Fi/Et: CO ₂ , O ₂ , N ₂ O, 5AA. RR, MAC Selection of waveform: CO ₂ , O ₂ , anesthetic agents


**Expert technologies:**

- Low sample flow capnography in the sidestream.
- Neonatal Oxygen-Cardio-Respirogramm.
- Pulse oximetry Masimo SET®.
- Anesthetic agent module (multigas).

Basic functions

Touchscreen	12.1" or 15"
Power supply	220 V/50 Hz
Battery	2 h of operation
Trends	digital and graphic, 240 h
Interface	wired or wireless connection to the central monitoring station (up to 32 monitors) wired connection to the hospital network USB for transferring patient data to PC
Alarm system	three levels of priority setting alarm limits alarm history
Data printing	three channels thermal printer
Medical calculators	for drug doses for oxygenation and ventilation parameters for renal function for hemodynamics parameters

Benefits

Low sample flow capnography in the sidestream	Necessary for neonatal patients. 50 ml/min sample volume is because of a small measuring chamber. Low sample flow improves capnography reliability in case of low respiratory flow (newborn infants and children)
Neonatal Oxygen-Cardio-Respirogramm	This is the simultaneous recording of heart rate, respiration rate, the blood oxygen saturation from newborn infants. Oxygen-Cardio-Respirogramm shows respiratory function, and maturity of the brain
Pulse oximetry Masimo SET®	The use of Masimo SET® pulse oximetry allows to assess real value of SpO ₂ in cases of disturbed perfusion or high levels motor activity of the patient
	Also, Masimo SET® pulse oximetry provides data to estimate the perfusion index. It is important in assessment of capillary blood flow

Treaton

www.treat-on.com

e-mail: info@treat-on.com

tel.: +7 (343) 304-60-57, +7 (343) 304-60-53

620100 Ekaterinburg, Russia, Sibirsky Tract 12/5

Triton Electronic Systems Ltd.



Treaton continuous innovation

Product Catalog

- Intensive Care Ventilators
- High Frequency Jet Ventilator
- Patient Monitors
- Depth of Anesthesia and Sedation Monitor